

US Heavy Duty Washer

M223/4-5

Selected and certified for AATCC Test methods requiring repeated home laundering (methods 88 B+C, 124, 130, 135, 142, 143, 150, 172 and 179). Available in US domestic and international electrical specifications.

SIZE	
WEIGHT	75 kg / 165 lbs
DIMENSIONS	699 mm (w) x 686 mm (d) x 1067 mm (h) / 28 in (w) x 27 in (d) x 42 in (h)

ORDERING INFO	
M223/4	US Heavy Duty Washer 120V / 60 Hz
M223/5	US Heavy Duty Washer 220 V / 50 Hz

STANDARDS
AATCC 88B+C, 124, 130, 135, 142, 143, 150, 159, 172, 179, 188

